

◎ **Over the Wire Design**


- Enhanced Pushability
- Outstanding Tractability
- Minimized Crossing Profile

◎ **Soft low profile tip**


◎ **Semi-Compliant Balloon**

- High Pressure
- Precision Placement
- Rapid Inf./Deflation


ValvuloPro® at aortic valve before expanding




ValvuloPro® expanded at aortic valve (3ATM)



Compliance Chart



Soft & low profile tip of ValvuloPro® vs. Competitor's at same balloon size(26mm)



Distributed in the US by: TT Medical USA, LLC.

📍 Duxbury MA 02332 ✉ sales@ttmedusa.com
 🌐 www.ttmedicalusa.com ☎ 508-269-2074

Manufactured by: Dongguan TT Medical, Inc.

Guangdong, PR China
 www.ttmedicalinc.com

Balloon Diameter (mm±10%)	Balloon Length(mm, ±10%)			Catheter Outer Dimension (mm)	Compatible Introducer Sheath (Fr)	Compliant Guidewire (inch/ mm)
	20	40	60			
4	VBC0420	VBC0440	VBC0460	1.32	4	0.018/ 0.46
5	VBC0520	VBC0540	VBC0560	1.50	5	
6	VBC0620	VBC0640	VBC0660			
7	VBC0720	VBC0740	VBC0760			
8	VBC0820	VBC0840	VBC0860			
9	VBC0920	VBC0940	VBC0960	2.00	7	
10	VBC1020	VBC1040	VBC1060			
11	VBC1120	VBC1140	VBC1160			
12	VBC1220	VBC1240	VBC1260			
13	VBC1320	VBC1340	VBC1360	2.50	8	
14	VBC1420	VBC1440	VBC1460			
15	VBC1520	VBC1540	VBC1560			
16	VBC1620	VBC1640	VBC1660			
17	VBC1720	VBC1740	VBC1760			
18	VBC1820	VBC1840	VBC1860	2.65	9	
19	VBC1920	VBC1940	VBC1960			
20	VBC2020	VBC2040	VBC2060			
21	VBC2120	VBC2140	VBC2160			
22	VBC2220	VBC2240	VBC2260	3.00	10	
23	VBC2320	VBC2340	VBC2360			
24	VBC2420	VBC2440	VBC2460			
25	VBC2520	VBC2540	VBC2560			
26	VBC2620	VBC2640	VBC2660			
27	VBC2720	VBC2740	VBC2760			
28	VBC2820	VBC2840	VBC2860			
29	VBC2920	VBC2940	VBC2960	12	14	
					15	